


VERSION FOR SIEMENS DIMENSION VISTA I500 AUTOMATED ANALYZER




DESCRIPTION:

Xylossay[®] is a reliable and automatable method for the quantification of xylose present in urine. **Xylossay**[®] is the recommended method for xylose detection after the administration of **LacTEST 0.45 g**.

METHOD RATIONALE:

The enzyme xylose dehydrogenase (XylDH)¹ catalyzes the specific oxidation of xylose to yield xylonolactone². As XylDH is a NAD⁺ dependent enzyme, this reaction requires a concomitant reduction of this cofactor to NADH:



Production of NADH can be spectrophotometrically detected and quantified by the increase of the absorbance at 340 nm during the reaction. This way, absorbance increments will be directly proportional to the amount of xylose present in the sample.

COMPONENTS OF THE KIT:

Xylossay[®] is supplied in a format suitable for performing 110 xylose determinations. Each kit contains a vial of each one of the following components:

VIAL	CONTENT AND FORMAT	AMOUNT
Brown	Phosphate buffer, 50 mM, pH 8.0 (Solution)	22 mL
Orange	β-Nicotinamide Adenine Dinucleotide (NAD ⁺) (Lyophilized)	66.34 mg ± 12%
Purple	D-Xylose Dehydrogenase (Lyophilized)	0.6 mg ± 10%
Pink	Calibrator: D-(+)-Xylose (Solution)	2 mL; 3.75 mg/dL ± 5%.

1. Sánchez-Moreno, I. et al. *J. Biotechnol.* (2016) 234:50-57.

2. Stephens, C. et al. *J. Bacteriol.* (2007) 189:2181-2185.

NON SUPPLIED ADDITIONAL EQUIPMENT:

Preparation of the reagents requires the use of micropipettes or systems suitable for dispensing the specified volumes. Vortex agitators can be used when indicated. Xylose quantification controls.

PREPARATION OF REAGENTS:

IMPORTANT: the preparation of the reagents must be conducted in the following order:

1. Dissolve the content of the purple vial (Xylose dehydrogenase, lyophilized) in **2.5 mL** of buffer (Phosphate buffer 50 mM, pH 8.0; brown vial). Mix gently to avoid loss of activity in the suspended enzyme. Keep cold during use if possible. This solution will be called **REACTIVE 2**.

2. Add 2 mL of buffer (Phosphate buffer 50 mM, pH 8.0; brown vial) in the vial that contains the lyophilized NAD⁺ (orange vial) and mix vigorously until its complete dissolution (vortex agitators can be used).

3. Add the totality of the 2 mL of dissolved NAD⁺ (orange vial + buffer) in the buffer vial (brown vial). Thus, the brown vial will contain **19.5 mL** of phosphate buffer plus NAD⁺ at the appropriated concentration. Keep cold if possible. This solution will be called **REACTIVE 1**.

4. The solution of xylose (Calibrator: standard xylose solution, pink vial) is ready for use.

CONSERVATION AND STABILITY:

Xylossay[®] must be stored at temperatures between 4° and 8°C until use, being stable under these conditions for at least 18 months. After being prepared, **REACTIVES 1** and **2** are stable for at least 39 days when maintained at temperatures between 4 and 10°C (stability in automated equipments). Both **REACTIVE 1** and **REACTIVE 2** can be frozen at -20°C without any activity loss. Both reagents can be frozen and thawed at least 6 times without altering their function.

IMPORTANT: to use frozen reactivities, these must be completely thawed (at room temperature). Once completely thawed, they must be agitated for total homogenization before use.

ASSAY PROTOCOL:

SAMPLE PREPARATION: Collection of urine samples must be performed strictly following the specifications indicated in **LacTEST 0.45 g** summary of product characteristics. Centrifugation of urine samples is recommended (according to the protocol established in each laboratory) to remove any possible precipitate that could affect the assay. Frozen samples must be agitated after thawing for their homogenization.

PARAMETERS: wavelength 340 nm, using the standard device temperature (usually 37°C). The instrumental application for testing must be designed according to the following specifications.

CALIBRATION AND CONTROL PROCEDURES: Calibration of xylose determination in the analyzers (**CAL**) must be carried out using the standard xylose solution supplied with the kit (Calibrator, pink vial) as concentration reference (3.75 mg/dL). For calibration and xylose quantification control it is recommended the use of at least two certified and validated levels of control: a low level, with xylose concentration between 0.8 and 3 mg/dL, and a high level,

with xylose concentration between 5 and 8 mg/dL. These controls are not supplied with the kit, and must be purchased independently.

STANDARD ASSAY PROTOCOL:

The shown volumes can be increased or reduced if the final proportion between them is maintained.

Reactive	Blank	Sample	CAL
Distilled water	24 µL	-	-
REACTIVE 1	100 µL	100 µL	100 µL
Sample	-	24 µL	-
Calibrator	-	-	24 µL
Mix the reaction and incubate for 5 min. (A1)			
REACTIVE 2	14 µL	14 µL	14 µL
Mix the reaction and incubate for 5 min. (A2)			

FINAL NUMBER OF MEASUREMENTS WITH EACH KIT:


Due to the dead volume of some automated analyzers, the number of determinations can be reduced. The next formulas allow calculating the final number of assays in function of the dead volume:

$$REACTIVE\ 1 = \frac{(19.5\text{ mL} - \text{Dead volume in mL}) \times 1000}{\text{Volume of REACTIVE 1 in the assay } (\mu\text{L})}$$

$$REACTIVE\ 2 = \frac{(2.5\text{ mL} - \text{Dead volume in mL}) \times 1000}{\text{Volume of REACTIVE 2 in the assay } (\mu\text{L})}$$

CALCULATIONS:

Absorbance values (340 nm) per assay:



A1 = Initial absorbance of the mixture REACTIVE 1 + Sample (incubated for 5 min).

A2 = Final absorbance after adding REACTIVE 2 (incubated for 5 more minutes, total time of reaction: 10 min).

Differences between the two absorbance values will be proportional to xylose concentration, which can be calculated using the standard xylose solution provided with the kit (Xylose standard solution, pink vial):

$$\Delta\text{Absorbance (340 nm)} = \Delta\text{Abs} = A2 - A1$$

$$\text{Sample concentration} = [\text{Sample}] \text{ (mg/dL)}$$

$$\text{Xylose concentration (Calibrator)} = 3.75 \text{ mg/dL}$$

$$[\text{Sample}] = \frac{\Delta\text{Abs (Sample)}}{\Delta\text{Abs (Calibrator)}} \times 3.75 \text{ mg/dL}$$

The total amount of xylose in the Sample (mg) can be calculated using the total volume of urine collected during the test.

$$\text{Xylose (mg)} = [\text{Sample}] \times \text{Volume of Sample (dL)}$$

EXAMPLE:

Determination of the total amount of xylose in the urine sample of a patient:

Volume of Sample (urine) = 557 mL = 5.57 dL

- **Sample assay:** A1 = 0.122

A2 = 0.146

$$\Delta\text{Abs (Sample)} = 0.146 - 0.122 = 0.024$$

- **Calibrator assay:** A1 = 0.0870

A2 = 0.166

$$\Delta\text{Abs (Calibrator)} = 0.166 - 0.0870 = 0.079$$

$$[\text{Sample}] = \frac{0.024}{0.079} \times 3.75 \text{ mg/dL} = 1.139 \text{ mg/dL}$$

$$\text{Xylose (mg)} = 1.139 \frac{\text{mg}}{\text{dL}} \times 5.57 \text{ dL} = 6.34 \text{ mg}$$

REFERENCE VALUES

Normal values in adults:

Test	Xylose in urine (mg)
LacTEST 0.45 g	≥ 19.18

Values below 19.18 mg indicate hypolactasia.

ANALYTICAL PARAMETERS:

The analytical parameters of **Xylossay**® application were determined during the validation of the technique in three automated analyzers (Roche Cobas c502, Werfen ILab 600 and Siemens Dimension Vista I500 analyzers).

Linearity: until at least 15 mg/dL in all the cases.

Limit of detection: between 0.13 and 0.49 mg/dL.

Reproducibility: The coefficients of variation between samples cannot be higher than 15%.

Accuracy: ≥ 88 %.

Carry-over: ≤ 4 %

ANALYTICAL PERFORMANCE:

The analytical parameters established for **Xylossay**® can slightly vary in function of the characteristics of the equipment used. Presence of some sugars as L-arabinose or glucose at high concentrations can affect xylose determination. For this reason, sample collection requires a fasting period of 10 hours. Under these conditions, no interference has been found in the collected urine.



Immunostep S.L.

Avda. Universidad de Coimbra, s/n
 Cancer Research Center (CIC)
 Campus Miguel de Unamuno
 37007 Salamanca (Spain)
 Tel. (+34) 923 294 827
www.immunostep.com

ANNEX – APPLICATION PARAMETERS

The following parameters were optimized in a Siemens Dimension Vista 1500 automated analyzer:


Advanced menu > Configuration > User-Defined Methods > Create New

Method name ID Units Mode

Execution | Reagent | Calculation | Calibration | Sample

Delivery	Time	Component 1	Remix	Component 2	Remix	Chase	Total Volume	Mix
D1	-21	R1 100 μL	None	--- 0 μL	None	0 μL	100 μL	None
S1	0	S 24.0 μL				0 μL	24 μL	Moderate
D2	392	R2 14 μL	None	--- 0 μL	None	0 μL	14 μL	Moderate
D3	39	--- 0 μL	None	--- 0 μL	None	0 μL	0 μL	None
S2	662	S 0.0 μL				0 μL	0 μL	None

Total Volume 138 μL



Reagent cartridge for User-defined tests in Dimension Vista 1500 equipment from SIEMENS (SIEMENS EMPTY Flex®, SIEMENS reference: KS999). In the case of XYLOSSAY®, Reagent 1 must be loaded in positions 1-5, and Reagent 2 in position 12. To add the reagents, the film that covers the wells of each position must be perforated in two of its corners, as indicated in the figure insert.

	Reagent	Tests	Life (hours)	Volume	
1	R1	27	2160	2992	Onboard Life (hours) 2160
2	R1	27	2160	2992	
3	R1	27	2160	2992	
4	R1	27	2160	2992	
5	R1	27	2160	2992	
6	Empty	0	0	0	
7	Empty	0	0	0	
8	Empty	0	0	0	
9	Empty	0	0	0	
10	Empty	0	0	0	
11	Empty	0	0	0	
12	R2	131	2160	2186	

Format
 Absorbance

Measuring
 340

Blanking

Time
 P1 310
 P2 620
 P3 -21
 P4 -21

Dilution
 1.00

IOD --- 0
FOD --- 0
FOD --- 0
FOD --- 0

Name	Formula
mau1	R93[340]
mau2	(r179[340])
CalculationNumber	mau2-mau1

Errors

Check --- **Read 1** -21 - **Read 2** -21 **Trigger** --- 0

Absorbance Add

Std. Curve **LINEAR** Interval (days) **30**

Calibrator Levels

2

Level	Weight	Replicates
1	1.00	3
2	1.00	3

Calibrator Dilutions

Create	Component	Diluent	Dilution Factor	Minimum Aspiration
<input type="checkbox"/>	L1			200
<input type="checkbox"/>	L2			200

Urine

Assay Range **Low** 0.2 **High** 15

Reference Range 0 15

Expiration (minutes) 120

Predilute

Sample volume 0

Dilute Factor 0

Diluent **SDIL**

Autodilution

Above Below

Sample volume 15 70

Dilute Factor 5 2

Diluent **SDIL**

Advanced menu > Calibration > Calibrators > New

Product Name	XILOCal	Bottle Values			Calibration Point Values				
Open Channel	YES		A	B	Unit	Level	1	2	Unit
Lot Number	XXXXX	Volume	2000	2000	µL				
Expiration Date	XXXX-XX-XX	XILO	0	3.75	mg/dL	XILO	0	3.75	mg/dL
OnBoard Stability	720								
Open Vial Stability	720								
Max Punctures	720								
Fluid Type	Urine								
Obsolete	NO								
		<input type="checkbox"/> Calibrator prepared externally							